Adm	Name	Class
Aum	! \ume !	C <i>uss</i>

GATITU SECONDARY SCHOOL P.O.BOX 327-0130 GATUNDU FORM THREE CCHEMISTRYPAPER 3 TERM 3 2014.

Question number	Maximum score	Students score
1	16	
2	12	
3	10	
TOTAL SCORE	40	

You are provided with

- (a) Solution M which is 2M sodium hydroxide.
 - (b) Solution N which is a dibasic acid whose concentration is 9.8gdm^{-3} with a general formula H_2A .
 - (c) Apparatus for titration
 - (d) Distilled water

You are required to

- I. Prepare solution P
- II. Determine molar mass of the dibasic acid in solution N.
- III. Determine formula mass of the acid radical (A) in the dibasic acid in solution N.

Procedure I

Using a pipette, transfer 25cm³ of solution M into the 250cm³ volumetric flask. Add distilled water to make up the volume to 250cm³. Swirl the volumetric flask five times to make a uniform solution. Label this solution as solution P. Use it for procedure II below.

(a) Calculate the concentration of solution P, (which you have prepared) in moles per dm³ (i.e molarity of solution P). (2mks)

Procedure II

Fill the burette with solution N.

Pipette 25 cm³ of solution P into a 250ml conical flask. Add three drops of phenolphthalein indicator. Titrate with solution N from the burette. Record your final and initial burette reading in the table below and complete the table. Repeat the titration two more times.

Titration	I	II	III
Final burette reading (cm ³)			
T-22-11 11 11 11 11 11 11 11 11 11 11 11 11			
Initial burette reading (cm ³)			
Volume of solution A used			
(cm ³)			

(a) Calculate the average volume of solution N used.

(4mks) (1mk)

- (c) Calculate
- (i) the number of moles of sodium hydroxide in the 25cm³ of solution P.

(2mk)

(ii) the number of moles of the dibasic acid two moles of sodium hydroxide reacts with		used given that (1mk)
(iii) the number of moles of the dibasic ac	oid in 1000cm ³ of solution N .	(2mk)
(iv) the relative formula mass of the dibasic (NB. concentration of the dibasic acid was		(2mk)
(v) the formula mass of the acid radical A	in the formula H_2A .	(2mk)
Question 2. Qualitative analysis (inorty) You are provided with solid Q. Carry observations and inferences in the space	out the experiment in tables below.	Write your
(a) Note the appearance of solid Q.		
Observation	Inferences	
(1mk)	(1mk)	

(b) Put all the solid Q in a clean dry test tube. Add about 10cm ³ the test- tube well. Retain this solution for the next tests.	of distilled water. Sl	hake

Observation	Inferences	
(1mk)	(1mk)	

- I. Divide the resulting solution into three portions.
- (i) To the first portion, add 2M sodium hydroxide dropwise till in excess.

Observation	vi sodium hydroxide dropwise till in exces	J.
	Inferences	
(11-)		
(1mk)	(1mk)	

(ii) To the second portion add 2M ammonia solution dropwise till in excess.

(iii) To the third portion, add three drops of 0.5M Barium nitrate solution. Record your observations. Add excess 2M nitric acid to the mixture and shake.

Observations Inferences			
Observations	Inferences		
	•		
	•		
(1 mls)			
(1mk)	(1mk)		
	(Time)		

(iv) To the fourth portion add three drops of lead nitrate solution. Add to the mixture about 5cm³ of distilled water shake and warm the mixture.

Observations	Inferences
(1mk)	(1mk)

(v) To the fifth portion add three drops of potassium iodide solution

Observation	Inferences	
	·	
/4 4\		
(1mk)	(1mk)	

- II. Write the formula of the cation and anion present in solid Q

Cation	(1mk)
Anion	(1mk)

Question 3. Qualitative analysis (Organic compounds)

You are provided with solid R. Carry out the experiment in tables below. Write your observations and inferences in the spaces provided.

(a) Put about one—third of solid R on the end of a clean metallic spatula and heat it with a Bunsen burner's flame.

	Inferences
servation	
	(11)
(1mk)	(1mk)
· · · · · · · · · · · · · · · · · · ·	f solid R in clean boiling tube. Add about 10cm ³ of ivide the solution into four portions . rops of bromine water
	Inferences
Observations	Interview
-	
	1
(1mk)	(1mk)
(i) To the second portion add	(1mk) d 3 drops of acidified potassium manganate (VII). Inferences
—	d 3 drops of acidified potassium manganate (VII).
(i) To the second portion add	d 3 drops of acidified potassium manganate (VII).
(i) To the second portion add	d 3 drops of acidified potassium manganate (VII).
(i) To the second portion add	d 3 drops of acidified potassium manganate (VII).
(i) To the second portion add	d 3 drops of acidified potassium manganate (VII). Inferences
(i) To the second portion add	d 3 drops of acidified potassium manganate (VII).
(i) To the second portion add Observations	d 3 drops of acidified potassium manganate (VII). Inferences
(i) To the second portion add Observations (1mk) (ii) To the third portion add	d 3 drops of acidified potassium manganate (VII). Inferences (1mk)
(i) To the second portion add Observations (1mk)	d 3 drops of acidified potassium manganate (VII). Inferences (1mk) d a spatula end full of sodium hydrogen carbonate
(i) To the second portion add Observations (1mk) (ii) To the third portion add	d 3 drops of acidified potassium manganate (VII). Inferences (1mk) d a spatula end full of sodium hydrogen carbonate
(i) To the second portion add Observations (1mk) (ii) To the third portion add	d 3 drops of acidified potassium manganate (VII). Inferences (1mk) d a spatula end full of sodium hydrogen carbonate
(i) To the second portion add Observations (1mk) (ii) To the third portion add	d 3 drops of acidified potassium manganate (VII). Inferences (1mk) d a spatula end full of sodium hydrogen carbonate

(iii) Dip a blue and red litmus papers in the fourth portion

Observations	Inferences
(1mk)	(1mk)