



231/3

BIOLOGY

PAPER 3

(PRACTICAL)

CONFIDENTIAL

Each candidate will require the following materials, reagents and equipment:-

1. 3 test tubes.
2. Hot water bath/means of heating.
3. Test tube rack.
4. 10ml measuring cylinder.
5. 4 droppers.
6. Test tube holder.
7. 10ml suspension X.
8. Filter paper.
9. Distilled water.
10. 1% copper (II) sulphate solution.
11. 10% sodium hydroxide.
12. Benedict's solution.

To prepare X:

Extract egg albumen from five eggs into a beaker. Add 20ml of olive oil then add distilled water to make 500ml suspension.