

3.7.3 Chemistry Paper 3 (233/3)

1. You are provided with:

- **Solution A**, 0.5 M copper(II) sulphate
- **Solid B₁**, metal B₁ powder
- **Solid B₂**, Iron powder
- **Solution C**, 0.02 M acidified potassium manganate(VII)

You are required to determine the:

- Enthalpy change for the displacement reaction between metal B₁ and copper(II) sulphate.
- Mass of iron that reacts with copper(II) sulphate in the displacement reaction.

PROCEDURE I

- (a) (i) Using a pipette and a pipette filler, place 25.0 cm³ of **solution A** into a 100 ml plastic beaker. Allow to stand for about 1 minute and then measure the temperature of the solution. Record the reading in **Table 1** as the initial temperature. Add all of **solid B₁** to the solution. Stir the mixture carefully with the thermometer and measure the highest temperature reached. This will take about **5 minutes**. Record the reading in **Table 1** as maximum temperature reached.

Table 1

Maximum temperature reached (°C)	
Initial temperature (°C)	
Change in temperature, ΔT_1 (°C)	

(3 marks)

(ii) Calculate the:

I number of moles of copper(II) sulphate used.

(1 mark)

II enthalpy change for the reaction of metal B₁ with one mole of copper(II) sulphate.

(Assume that for the mixture, specific heat capacity = 4.2 Jg⁻¹K⁻¹ and density = 1.0 g cm⁻³)

(1 mark)



- (b) Repeat **procedure I, (a) (i)** with all of **metal B₂** (iron powder) in place of **metal B₁**. The maximum temperature is reached after about 8 minutes. Record the temperature readings in **Table 2**. **Retain the mixture for use in PROCEDURE II.**

Table 2

Maximum temperature reached (°C)	
Initial temperature (°C)	
Change in temperature, ΔT_2 (°C)	

(3 marks)

- (c) Compare the changes in temperature ΔT_1 and ΔT_2 and comment on the differences.

(2 marks)

PROCEDURE II

- (i) Fill a burette with **solution C**.
- (ii) Filter the mixture obtained in **procedure I (b)** into a 250 ml volumetric flask. Wash the residue with distilled water and add into the flask. Add more distilled water to make up to the mark. Label this as **solution B₂**.
- (iii) Using a pipette and a pipette filler, place 25.0 cm³ of **solution B₂** into a 250 ml conical flask. Titrate **solution B₂** with **solution C** until a permanent pink colour just appears. Record the readings in **Table 3**.


Repeat step (iii) and complete **Table 3**.

- (d) **Table 3**

	I	II	III
Final burette reading			
Initial burette reading			
Volume of solution C used, cm ³			

(4 marks)

- (e) Calculate the average volume of **solution C** used. (1 mark)
- (f) The equation for the reaction between manganate(VII) and iron(II) ions is:



Calculate the number of moles of:

- (i) potassium manganate(VII) used. (1 mark)



- (ii) iron (II) ions in 25.0 cm³ **solution B** (1 mark)
- (iii) iron that reacted with copper(II) sulphate. (1 mark)
- (g) Determine the mass of iron that reacted. (RAM of Fe = 55.8) (1 mark)

2. You are provided with:

- **Solid K**
- Aqueous ammonia
- Aqueous sodium sulphate
- Dilute nitric(V) acid
- Wooden splint

Solid K is suspected to be **lead(II) carbonate**.

- (a) From the reagents provided, select and describe **three** tests that could be carried out **consecutively** to confirm if **solid K** is lead(II) carbonate. Write the tests and expected observations in the places provided.

(i)

Test 1	Expected Observations

(1 mark) (1 mark)

(ii)

Test 2	Expected Observations

(1 mark) (1 mark)

(iii)

Test 3	Expected Observations

(1 mark) (1 mark)

- (b) Carry out the tests described in (a) using **solid K** and record the observations and inferences in the spaces provided.

(i) Test 1

Observations	Inferences

(½ mark)

(½ mark)

(ii) Test 2

Observations	Inferences

(1 mark)

(2 marks)

(iii) Test 3

Observations	Inferences

(1 mark)

(1 mark)

3. You are provided with an organic compound **solid M**. Carry out the following tests. Record the observations and inferences in the spaces provided.

- (a) Place **all** of **solid M** in a boiling tube. Add about 10 cm³ of distilled water and shake. Retain the solution for use in procedure (b) (i), (ii) and (iii).

Observations	Inferences

(1 mark) (1 mark)

- (b) Use about 2 cm³ portions of the mixture in a test tube for tests (i), (ii) and (iii).

- (i) To the first portion, add all the solid sodium carbonate provided.

Observations	Inferences

(1 mark) (1 mark)

- (ii) To the second portion, add two drops of acidified potassium manganate(VII) and warm the mixture.

Observations	Inferences

(1 mark) (2 marks)

- (iii) To the third portion, add about 2 cm³ of acidified potassium dichromate(VI). Heat the mixture to boiling and allow to stand for about 2 minutes.

Observations	Inferences

(1 mark) (1 mark)

